



**PRISMA**  
COLLECTION

# BESPOKE PROGRESSIVES

## Advanced Lens Design

An Optometrist's Guide to Bespoke Multifocal,  
Occupational & Specialised Lens Designs



SWISS  
SMART  
INNOVATION

A Product of  
**VISION8**  
GLOBAL



# FIND US AT THESE **LOCATIONS**



## **CAPE TOWN**

27 Moody Ave,  
Epping 1,  
7475



## **JOHANNESBURG**

9 Delphi St,  
Eastgate, Sandton,  
2090





## **DURBAN**




22 Hillclimb Rd,  
Westmead,  
3610

# TABLE OF CONTENTS

## General Information

 DAFD - Dynamic Adaptable Frame Design	Pg 4-7
 DVS - Dynamic Vision System	Pg 8-9
 BES - Balance Eye Smoothing	Pg 10-11




## Progressive Designs & Centration Charts






 Essential Q Multifocal Starter Level	Pg 12-13
 InfinityFLEX Multifocal Premium	Pg 14-15
 Magic Prime 8k Multifocal Elite	Pg 16-17

## Specialised Progressive Design & Centration Chart

 Max Middle	Pg 18-19
--	----------

## Specialised Occupational Designs & Centration Charts






 Office Smart	Pg 20-21
 DRIVEmax	Pg 22-23
 Sports Agility	Pg 24-25

 EEZ Anti-fatigue	Pg 26-27
 Bifocal Blend 4	Pg 28-29
 TIB28 True Invisible Bifocal	Pg 30-31
 Centration Charts - Bifocal Blend 4 & True Invisible Bifocal	Pg 32-33
 MyopiaX	Pg 34-37


## Lens Design Summary Chart

Pg 38-39

## Other Lab Services

 Proctoguard Coatings	Pg 40-43
 Tintech Lens Tinting	Pg 44-45
 Eqlipsar Photochromatic	Pg 46-49
 Lens Thickness Comparison	Pg 50
 Products, Services & SAOA Codes	Pg 51






# BESPOKE PROGRESSIVES

The **Prisma Luxury Collection, by OTL Lens Lab**, is designed to cater to the unique visual needs of every individual, offering exceptional clarity, expansive fields of vision, and minimal distortion.

Each lens in the Prisma collection is **meticulously crafted to match the specific curvature of the wearer's eyes**, ensuring a perfect fit and optimal visual performance. Whether for daily activities, sports, or specialized occupational tasks, the Prisma Luxury Collection provides unparalleled visual clarity and comfort.



SWISS  
SMART  
INNOVATION



Experience the forefront of innovation with our US patented DAFD Technology  
- the state-of-the-art in freeform lens design

## Dynamic Adaptable Frame Design

Due to geometrical aspects, not all modern and fashionable frame shapes can be used with progressive designs. Although frame shapes are a starting point for buying a new pair of eyeglasses, they are likely to become a limitation when choosing the progressive designs.

In most cases, non-conventional frame shapes prevent progressive zones from being properly defined, and vision comfort is severely affected. In such cases the customer has to change the selected frame shape with a more suitable one. The main challenge today is to accomplish customer fashion frame requests with the best optical vision performance possible.




Dynamically sets the far and near zone length according to the frame shape

**Adaptable far and near width**

Designed to move with the frame shape, ensuring a flawless fit

**Dynamic decentration**



**Frame design rotation**

Aligns progressive zones in every design to fully utilise the frame's visual area

# Prisma DAFD: The Self-Designing Progressive Lens with AI Technology

The foundation of Prisma DAFD (Dynamic Adaptable Frame Design) lies in Prisma's DVS (Dynamic Vision System) Technology, enabling the individual customization of design settings, including far and near zones, inset, decentration, and more. By combining DVS and Prisma AI technology, along with the powerful Prisma Cloud computing infrastructure, DAFD design is born. Frame shape data sent to Prisma LDS automatically replicates the exact frame shape in the cloud and adjusts the progressive design to match the selected frame.

- **Automatic Design Rotation**

Aligns progressive zones to fully utilise the frame's visual area

- **Automatic Design Decentration**

Guarantees optimal placement for visual clarity


- **Dynamic Progressive Zone Adjustment**

Automatically customizes the width of far and near vision zones (DVS Technology)

- **AI-Driven Optimization**


Powered by Prisma AI, the system analyzes all possible configurations, cross-referencing historical customer data to select the most effective design

This ensures the perfect balance between reduced astigmatism, wider clear vision areas, and optimal lens thickness.






An advanced self-designing progressive lens engineered to adapt to any frame shape, creating a personalised visual solution that aligns with the wearer's unique vision requirements and style preferences.




### Traditional Prisma Progressive Design



Manual design settings

### Prisma Progressive Design with DAFD Technology



Automatic near vision rotation

Optimized peripheral areas

Expanded far vision area



## DVS – Dynamic Vision System

**One of Prisma's pioneering technologies, serving as the foundation for all other solutions introduced to the market.**

The DVS (Dynamic Vision System) Technology was developed to address one of the key challenges in the optical industry: the customization of progressive lenses. Currently, there are numerous specific designs aimed at meeting various user needs, such as prolonged outdoor use, intensive reading, or frequent use of digital devices. However, Prisma has revolutionized the traditional approach by creating a single progressive design capable of dynamically and precisely adapting to the individual needs of every wearer. By leveraging advanced mathematical algorithms, the DVS Technology enables adjustments to the width of the main vision areas – distance, near, and intermediate – within the same progressive design.

- Extension of the near zone: ideal for users who frequently use smartphones or tablets.
- Extension of the far zone: designed for drivers or those who primarily engage in outdoor activities.
- Softness of the design: for a comfortable and natural visual experience.

Adjustments to the vision zones must also consider peripheral astigmatism, a critical factor to balance alongside the length of the progressive channel. Prisma has made it possible to integrate infinite real world combinations into a single progressive design, thanks to the DVS Technology applied to various mathematical models, such as those based on hyperbolas or spline tensors.







## Extreme Flexibility and Adaptability

The DVS Technology provides unprecedented flexibility, adapting even to specific or extreme requests and ensuring rapid adaptability for the wearer. This means using a single design to meet the diverse needs of optometrists. With the Dynamic Vision System, Prisma delivers a unique solution combining innovation, customization, and visual comfort; meeting the demands of a continuously evolving market.

“  
**Prisma has developed a patented technology by using mathematical algorithms that allow it to make changes to its width of the main areas of vision within its lens design. This intrinsic flexibility of the design offers the advantage of customizing the design by using the different width parameters based on the purpose of using the lens.**  
 ”

## Available Settings for Progressive Zones

			
<p><b>Far vision optimization</b></p> <p>Ideal for improving vision in the upper part of the lens.</p>	<p><b>Near vision optimization</b></p> <p>Designed for those who require a wider and more precise lower zone.</p>	<p><b>Wide frames</b></p> <p>Perfect balance between near, intermediate and far zones in a larger frame.</p>	<p><b>Small frames</b></p> <p>The values of near, intermediate and far move for optimal adaptation in a smaller or narrow frame.</p>



## BES – Balance Eye Smoothing

The problems related to the adaptation in the use of progressive lenses are mainly due to the aberrations present in both the progressive and peripheral areas of vision, These force the wearer to unnatural movements of the head.




" The Prisma designs, thanks to the mathematical methodologies applied natively in the design phase, have always guaranteed a sharpness of the progressive areas. (Far, near and channel) close to 99%. Perceived astigmatism is therefore tending to 0."


The poor sharpness limits of the peripheral vision areas are due to the integration of the different geometries of the distant and near areas, thus creating not only high astigmatism values, but also a lot of non-homogeneity of the cylinders. These fluctuations in values, when reading the maps, show themselves as "waves of the sea". In nature, the more waves there are in the sea, the more agitated it appears to us. The higher the waves, the more there is a storm. In progressive lenses it works the same way, the more sudden changes in astigmatism there are and the more you have vision difficulties, the higher the astigmatism values are and the more you experience lens wear problems, such as headaches and nausea.

The Prisma designs, thanks to the mathematical methodologies applied natively in the design phase, have always guaranteed a sharpness of the progressive areas (Far, near and channel) close to 99%. Perceived astigmatism is therefore tending to 0. In recent years, the Prisma innovation team has focused on optimizing vision in the peripheral areas and Balance Eye Smoothing (BES) Technology is the result of this research. The BES combines different mathematical techniques that allow you to keep the average values of power and the values of astigmatism under control, standardizing them as much as possible within the peripheral areas. In this way, average astigmatism values are obtained that do not disturb binocular vision.


The BES, now applied in all Prisma designs, has composed and harmonized the set of progressive surfaces, ensuring greater comfort and sharpness of vision in all directions, so as to reduce the adaptation time to the progressive lens and limiting - in some cases eliminating the disadvantages of aberrations in the peripheral areas.



Angular and uneven areas




Softening in the passages between near and far areas to avoid the tooth effect




### Traditional lenses

The different areas are presented in a disconnected and non-uniform way



### Prisma lenses

Thanks to our mathematical algorithms, progressive zones are balanced and smoothed



## Essential Q - Multifocal Entry Level

Discover the ideal solution for newcomers to progressive lenses who rely heavily on near vision: the Essential Q Design. Customizable to accommodate individual needs and frame styles, it delivers an unparalleled level of comfort and performance.



Unlock the power of personalized vision with Prisma DVS Technology (Dynamic Vision System), enabling precise adjustments to both far and near vision zones. Experience a seamless journey as the near vision area gracefully expands from the onset, ensuring effortless transitions between various distances.

Elevate your visual experience with the Essential Q Design, meticulously crafted to meet your individual visual needs and frame preferences. Experience unmatched precision and comfort, setting a groundbreaking standard within the basic line of lenses.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White



## Progressive Designs



**Essential Q**  
Multifocal Entry Level


FAR VISION ● ● ○

INT. VISION ● ● ○

NEAR VISION ● ● ○

DESIGN DETAILS

- CL 9mm  
↓  
17mm
- Inset →  
0-5mm
- ⊕ ⊖  
8  
MAX Diopter
- SOFT
- Far Vision  
40-100mm
- Reading Area  
10-40mm
- Min Fitting  
13-21mm




# Essential Q


Multifocal Entry Level

## Centration Chart



SWISS  
**SMART**  
INNOVATION





## Infinity FLEX - Multifocal Premium

For the dynamic individual seeking a dependable solution across various activities. Infinity FLEX provides clear central vision and exceptional peripheral vision. With versatile channels available, this design seamlessly adapts to diverse wearer preferences.


Infinity FLEX Premium Design, has undergone numerous updates over time, pushing the boundaries of design possibilities to the extreme.

With a maximum addition of up to 6 and an additional 2 extra long corridors, it opens up infinite possibilities for customization and usage combinations.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White

Engraving  
Symbol  
**FX**

## Progressive Designs



**Infinity FLEX**  
Multifocal Premium

**FAR VISION** DYNAMIC

**INT. VISION** ● ● ○

**NEAR VISION** DYNAMIC

**DESIGN DETAILS**

- CL 9mm  
↓  
17mm
- Inset → 0-5mm
- ⊕ ⊖  
8  
MAX Dioptr
- SOFT
- Far Vision 40-100mm
- Reading Area 10-40mm
- Min Fitting 13-21mm




**Infinity FLEX**  
Multifocal Premium

**Centration Chart**



SWISS  
**SMART**  
INNOVATION





## MagicPrime 8K - Multifocal Elite


The design originates from Prisma DVS (Dynamic Vision System) Technology, **enabling you to tailor the width of both far and near vision zones according to your specific requirements.** This innovation marks the debut of the Prisma Elite multipurpose progressive design.

In its standard setup, its minimal astigmatism ensures seamless compatibility with frames of substantial height.

Elevate your visual experience with our Elite progressive lens design, where versatility meets unparalleled performance. Choose excellence, choose customization, choose our dynamic vision solution for a truly tailored experience.

AVAILABLE INDEXES					Engraving Symbol <b>MP</b>
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White	

## Progressive Designs



**MagicPrime 8k**  
Multifocal Elite

**FAR VISION** DYNAMIC  
**INT. VISION** DYNAMIC  
**NEAR VISION** DYNAMIC

**DESIGN DETAILS**

- CL 9mm → 17mm
- Inset → 0-5mm
- MAX Dioptr 8
- SOFT
- Far Vision 40-100mm
- Reading Area 10-40mm
- Min Fitting 13-21mm




# MagicPrime 8k

Multifocal Elite

## Centration Chart



SWISS  
**SMART**  
INNOVATION





## MaxMiddle

Experience uncompromised vision: MaxMiddle delivers clarity across every angle. Discover the innovation behind seamless transitions for peripheral, intermediate, and central vision zones.


Standout clarity: explore the remarkable width of the intermediate zone, offering clear and comfortable vision at reading and medium to short distances. Engineered with cutting edge technology, this lens enhances visual perception for effortless focus.

The lens is crafted to minimize distortions and aberrations, ensuring a superior visual experience in all situations. The progressive correction of presbyopia is seamlessly integrated, allowing for a smooth transition between various vision zones.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White



## Specialised Progressive Design



**MaxMiddle**

- FAR VISION ● ● ● ◐
- INT. VISION ● ● ● ●
- NEAR VISION ● ● ● ◐

**DESIGN DETAILS**

- CL 11mm  
↓  
19mm
- Inset → 2.2-2.3mm
- Inset: ⊕ ⊖  
6  
MAX Diopter
- SOFT (feather icon)
- Far Vision: 78mm (horizontal double-headed arrow)
- Reading Area: 22mm (horizontal double-headed arrow)
- Min Fitting: 15-21mm (vertical double-headed arrow)




MaxMiddle

Centration  
Chart



SWISS  
SMART  
INNOVATION





## Office Smart


This design is purpose built for modern workstation demands, supporting sustained focus at intermediate distances while maintaining comfortable awareness of the surrounding environment. By carefully balancing the visual zones and softening the distance power, it reduces visual strain and enables effortless gaze shifts between screen, desk, and room - promoting comfort and visual efficiency throughout extended periods of use.

Regression Value	Ranges	Viewing Distance at Top
0.75D	Computer Range	Focus stays within 1 meter (Close desk work)
1.50D	Meeting Range	Focus up to 2 meters (Large monitors/Cubicle)
2.25D	Office Range	Focus up to 4 meters (Meeting rooms)

AVAILABLE INDEXES					Engraving Symbol	QR Code
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White		

SCAN QR CODE FOR DESKTOP DEPRESSION INFORMATION TABLE

## Specialised Occupational Designs



FAR VISION ● ○ ○  
 INT. VISION ● ● ●  
 NEAR VISION ● ● ●

DESIGN  
DETAILS

CL  
↓  
13mm

Inset  
→  
2mm

⊕ ⊖  
4  
MAX  
Diopter

SOFT

Intermediate  
Vision  
↔  
84mm

Reading  
Area  
↔  
24mm

Min  
Fitting  
18mm  
↑↓





**Office Smart**

**Centration Chart**



**SWISS  
SMART  
INNOVATION**





## DRIVEmax


Introducing DRIVEmax, a progressive lens design meticulously engineered for drivers and individuals desiring an extensive far vision area. Suitable for both low and high prescription additions. This design delivers unmatched clarity and comfort while driving.

Tailor your progressive vision experience with Prisma DVS Technology (Dynamic Vision System). This innovative feature enables personalized adjustments to the width of the near vision zone in Prisma driver's lenses, guaranteeing optimal visual performance customized to your specific requirements.

DRIVEmax features a 4° design rotation, enhancing the lateral view and providing improved visibility of the rear-view mirrors. This significant innovation benefits drivers without altering the technical parameters of the lens design.

AVAILABLE INDEXES					Engraving Symbol <b>DM</b>
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White	

## Specialised Occupational Designs



**DRIVEmax**

FAR VISION ● ● ●  
 INT. VISION ● ● ○  
 NEAR VISION DYNAMIC

DESIGN  
DETAILS

CL  
9mm  
↓  
17mm

Inset  
→  
0-5mm

⊕ ⊖  
8  
MAX  
Diopter

SOFT

HARD

Far  
Vision  
FULL

Reading  
Area  
10-40mm

Min  
Fitting  
13-21mm




**DRIVEmax**

## Centration Chart



SWISS  
**SMART**  
INNOVATION





## Sports Agility

Unlock unparalleled sports performance with Sports Agility progressive lenses, meticulously crafted for athletes pursuing unmatched visual clarity and precision. Engineered to exceed the rigorous demands of sports frames, these lenses provide unrivaled distance vision and sharp peripheral clarity, essential for achieving peak athletic excellence.


Powered by state-of-the-art Prisma technology, these lenses guarantee crystal clear vision throughout, elevating binocular vision for improved depth perception and precision. With an expansive far vision range, these lenses excel in intense activities, yet maintain clear sight of nearby surroundings as needed.

Tailored for wrap frames, Sports Agility lenses deliver uncompromising vision quality and comfort, enhancing every sporting pursuit.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White



## Specialised Occupational Designs



**Sports Agility**

- FAR VISION ● ● ●
- INT. VISION ● ● ◐
- NEAR VISION ● ● ○

DESIGN DETAILS

- CL  
↓  
13mm
- Inset  
→  
2.2mm
- ⊕ ⊖  
4  
MAX  
Diopter
- SOFT
- Far Vision  
↔  
90mm
- Reading Area  
↔  
24mm
- Min Fitting  
↕  
19mm




Sports Agility

## Centration Chart



SWISS  
SMART  
INNOVATION





## EEZ Anti Fatigue

Introducing our revolutionary EEZ Anti Fatigue Design, crafted specifically for students to combat eye strain and fatigue. Engineered for 24-hour wear, it utilizes cutting-edge technology to support and relax young eyes.


By mitigating the effects of screen time and study sessions, this design enables seamless transitions between tasks without causing strain.

Whether reading, studying, or using digital tools, the EEZ Anti Fatigue Design provides constant refreshment and revitalization, prioritizing well-being through clarity and comfort.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White



## Specialised Occupational Designs



- FAR VISION ● ● ●
- INT. VISION ● ● ◐
- NEAR VISION ● ● ●

DESIGN DETAILS

- CL: 13mm
- Inset: 2mm
- MAX Diopter: 4
- SOFT
- Far Vision: 120mm
- Reading Area: 24mm
- Min Fitting: 18mm




**EEZ  
Anti Fatigue**

**Centration  
Chart**



**SWISS  
SMART  
INNOVATION**





# Bifocal Blend 4



NEAR VISION WIDTH

## Generated on all blanks


**Free-form Bifocal Blend 4 lenses aim to eliminate this visible line, providing a seamless transition between the different vision zones.** This design not only enhances the aesthetics but also offers better visual comfort and clarity.

One of the key features of the Prisma Bifocal Blend 4 lenses is the ability to customize the near width of the lens according to the wearer's needs. Eyewear professionals can tailor the near width based on factors such as the individual's prescription, pupillary distance, reading habits, and other visual preferences.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White



# BIFOCALS



You can minimize the near vision area width to maximize the far vision area

or you can expand the near vision area width as you want

The minimum near area width (also know as "crescent") starts at 24mm diameter and the maximum width is set to use the same value of the lens diameter. Near area width personalization can be done with 1mm step. Near area vertical position (also referred as segment vertical position) is also customizable spanning from 4mm to 12mm.

CL  
↓  
8mm

Near Vision  
↔  
24mm to Lens Diameter width

SOFT

⊕ ⊖  
8  
MAX Diopter

Far Vision  
↔  
FULL



# TIB28 - True Invisible Bifocal

## Advanced Bifocal Lens with Transition Channel Technology

Traditional bifocal lenses present an abrupt optical discontinuity between the distance and near zones, often resulting in astigmatic jumps and wearer discomfort.

True Invisible Bifocal represents the next generation of bifocal lenses, featuring a **28mm near segment with precise optics and a smooth transition channel that eliminates visual jumps**. Its design allows for simplified polishing and a high quality finish, while providing comfortable, natural vision for all day use.



## LENS COMPARISON



**Traditional Bifocals:** High image jump, limited intermediate vision, moderate adaptation, visible segment line, medium polishing complexity.



**Progressive Lenses:** No image jump, peripheral astigmatism, no segment line, high polishing complexity.



**True Invisible Bifocal:** No image jump, minimal distortions, easy adaptation, fixed 28mm invisible segment, simplified polishing.

AVAILABLE INDEXES				
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White





Combines the appearance of a single-vision lens with the benefits and advanced technology of bifocal lenses

 **Distance Zone**  
Wide & distortion free

 **Segment (28mm)**  
Sharp, stable optics for near tasks

 **Transition Channel**  
Seamless shift between far & near



### OPTICAL CHARACTERISTICS

- No prism jumps during eye movement
- Minimal astigmatic distortions compared to progressive lenses
- Easy adaptation for all wearers



### MANUFACTURING ADVANTAGES

- Simplified surfacing and polishing process
- High reproducibility of results
- Consistently optimal surface finish



### CLINICAL & PRACTICAL APPLICATIONS

- Stable and distortion free transitions
- Suitable for wearers sensitive to progressive lens astigmatism
- Laboratories seeking efficient production with high optical quality

CL  
↓  
12mm

Near Vision  
  
28mm

Far Vision  
  
FULL


  
8  
MAX  
Diopter

*\*Kindly note:  
Measurements  
should be taken  
as a Multifocal*



# Bifocal Blend 4


## Centration Chart



NEAR VISION WIDTH



SWISS  
SMART  
INNOVATION






# TIB28 - True Invisible Bifocal

## Centration Chart




SWISS  
**SMART**  
INNOVATION





**MyopiaX**



## Clear vision with incredible details

The traditional Myopia lenses create a blurred image on the peripheral area



MyopiaX ensures clear vision across the entire lens area, allowing the focal point of incoming light to shift backwards so that it converges precisely on the retina through the centre of the lens



The circular defocus changes the focus points' positions while moving away from the lens centre to follow the retina eye curve, as accurately as possible. The result is a lens focusing the light on the retina in correspondence of every point of the lens without blurring the images.

Single vision design is used to control children myopia that combines effectiveness and portability. Prisma has developed MyopiaX Design, a next generation single vision freeform design with circular defocus limiting the worsening of myopia in young people.

MyopiaX is based on a clinically validated optical concept, supported by published scientific evidence and ongoing long-term clinical evaluation, combining efficacy with high wearability to support real-world myopia management in children. Clinical evidence shows greater effectiveness in younger children (6-9 years), supporting early intervention.

MyopiaX enables eye care professionals to:


- move from simple vision correction to active myopia management
- rely on a clinically validated optical concept
- offer a non-invasive, spectacle-based solution
- support real-world effectiveness through high wearability and compliance
- confidently recommend early intervention, especially in younger children because MyopiaX lenses are worn like standard glasses, they integrate naturally into everyday optometric practice.



## MyopiaX - Heat Map




SWISS  
**SMART**  
INNOVATION




**MyopiaX Toddler**  
5-8 years

MyopiaX Toddler with a total defocus of 2.5 Diopters for kids starting from 5 to 8 years old. Compared to a traditional single vision myopia lens, MyopiaX Toddler is 25% thinner.



**MyopiaX Junior**  
9-12 years


MyopiaX Junior with a total defocus of 1.6 Diopters for a youngster from 9 to 12 years old. Compared to a traditional single vision myopia lens, MyopiaX Junior is 18% thinner.



**MyopiaX Teen**  
13-17 years

MyopiaX Teen with a total defocus of 1.0 Diopters for teenagers from 13 to 17 years old. Compared to a traditional single vision myopia lens, MyopiaX Teen is 10% thinner.

## How it works - The MyopiaX Revolution



Power Map with axes directions

Myopia X has a fixed power central zone which diameter length is calculated from the real base value and the power value of the requested lens. Generally, diameter range from 9 to 12 mm. From the central zone to the lens border, MyopiaX has a power depression starting from 1 to 2.5 Dioptrics, depending on which design is used (MyopiaX Toddler, MyopiaX Junior, MyopiaX Teen).

Thanks to the ray tracing technique, MyopiaX design has a circular defocus and from the design center to the design peripheral area, power decreases and cylinders slowly increase. Cylinder axes are always directed to the lens center (radial direction) for better comfort on all eye rotations.

MyopiaX lenses can be produced on any material and in any diameter, with any blanks and any index on the market. Customers can specify any individual parameters like wrap, panto and back vertex distance. Sphere range from 0.00 to -20.00 and Cylinder max 6, depending on the age of the wearer, MyopiaX has 3 different design models and for each of

them thickness is significantly improved thanks to the circular parabolic defocus. MyopiaX surface is clear without porosity.

According to Clinical Evidence Published Outcomes (at 12 months): Compared to standard single vision lenses, the peripheral defocus design showed:

- ~30% reduction in myopia progression (SER)
- ~24% reduction in axial length growth
- Statistically significant results
- Greater effectiveness in younger children (6-9 years)
- High wearing compliance and visual comfort

These results demonstrate a clinically meaningful and safe myopia control effect.

AVAILABLE INDEXES					Engraving Symbol <b>M</b>
1.50 White	1.56 Sunglare	1.67 White	1.67 Sunglare	1.74 White	



# LENS DESIGN SUMMARY



SWISS  
SMART  
INNOVATION

LENS	ENGRAVING SYMBOL	TYPE	VISIONS	CORRIDOR LENGTH	MIN. FIT HEIGHT mm	MAX ADD	INSET	FAR VISION	READING AREA
<b>Progressive</b> EssentialQ Multifocal Entry Level	<b>EQ</b>	 SOFT	FAR VISION ● ● ○ INT. VISION ● ● ○ NEAR VISION ● ● ○	9mm ↓ 17mm	Min Fitting ↑ ↓ 13 to 21 mm	⊕ ⊖ 8 MAX Diopter	Inset ↔ 0 to 5	Far Vision  40-100mm	Reading Area  10-40mm
<b>Progressive</b> Infinity FLEX Multifocal Mid-Range	<b>FX</b>	 SOFT	FAR VISION DYNAMIC ● ● ○ INT. VISION ● ● ○ NEAR VISION DYNAMIC ● ● ○	9mm ↓ 17mm	Min Fitting ↑ ↓ 13 to 21 mm	⊕ ⊖ 8 MAX Diopter	Inset ↔ 0 to 5	Far Vision  40-100mm	Reading Area  10-40mm
<b>Progressive</b> MagicPrime 8k Multifocal Premium Range	<b>MP</b>	 SOFT	FAR VISION DYNAMIC ● ● ○ INT. VISION DYNAMIC ● ● ○ NEAR VISION DYNAMIC ● ● ○	9mm ↓ 17mm	Min Fitting ↑ ↓ 13 to 21 mm	⊕ ⊖ 8 MAX Diopter	Inset ↔ 0 to 5	Far Vision  40-100mm	Reading Area  10-40mm
<b>Specialised</b> MaxMiddle Wide Intermediate Range	<b>MD</b>	 SOFT	FAR VISION ● ● ○ INT. VISION ● ● ● NEAR VISION ● ● ○	11mm ↓ 19mm	Min Fitting ↑ ↓ 15 to 21 mm	⊕ ⊖ 6 MAX Diopter	Inset ↔ 2.2 to 2.3	Far Vision  78mm	Reading Area  22mm



# LENS DESIGN SUMMARY



SWISS  
SMART  
INNOVATION

LENS	ENGRAVING SYMBOL	TYPE	VISIONS	CORRIDOR LENGTH	MIN. FIT HEIGHT mm	MAX ADD	INSET	FAR VISION	READING AREA
Specialised Office Smart	OS	SOFT	FAR VISION ● ○ ○ INT. VISION ● ● ● NEAR VISION ● ● ●	↓ 13mm	Min Fitting ↑↓ 18mm	⊕ ⊖ 4 MAX Diopter	Inset ↔ 2mm	Far Vision ↔ 84mm	Reading Area ↔ 24mm
Specialised DRIVEmax	DM	SOFT HARD	FAR VISION ● ● ● INT. VISION ● ○ ○ NEAR VISION DYNAMIC	9mm ↓ 17mm	Min Fitting ↑↓ 13 to 21 mm	⊕ ⊖ 8 MAX Diopter	Inset ↔ 0 to 5	Far Vision ↔ FULL	Reading Area ↔ 10-40mm
Specialised Sports Agility	SA	SOFT	FAR VISION ● ● ● INT. VISION ● ● ● NEAR VISION ● ● ○	↓ 13mm	Min Fitting ↑↓ 19mm	⊕ ⊖ 4 MAX Diopter	Inset ↔ 2.2	Far Vision ↔ 90mm	Reading Area ↔ 24mm
Specialised EEZ Anti-fatigue	AF	SOFT	FAR VISION ● ● ● INT. VISION ● ● ● NEAR VISION ● ● ●	↓ 13mm	Min Fitting ↑↓ 18mm	⊕ ⊖ 4 MAX Diopter	Inset ↔ 2mm	Far Vision ↔ 120mm	Reading Area ↔ 24mm

# 7 PROTECTIVE FEATURES



**ANTI  
REFLECTION**



**SCRATCH  
RESISTANT**



**OIL & SMUDGE  
REPELLENT**



**ANTI-STATIC**



**WATER  
REPELLENT**



**DURABLE**




**ENHANCED  
DAY & NIGHT  
VISION**

T's & C's apply

**PROCTO™  
GUARD  
COATINGS**

**PROCTO™  
TANZANITE  
BLUE BLOCK**

**PROCTO™  
JADE GREEN**



**(OPTIONAL)  
OLEOPHOBIC WITH  
SUPER HYDROPHOBIC**

**(OPTIONAL)  
HYDROPHOBIC  
WITH ANTI-STATIC**

**HYDROPHOBIC**

**6 AR MULTI NANO-LAYER  
COATINGS**

**HARD COAT**

**SUBSTRATE (LENS SURFACE)**



**PREMIUM COATING FOR EXCEPTIONAL CLARITY & PROTECTION**



PROCTO™  
GUARD  
COATINGS



PRISMA  
COLLECTION



PROCTO™  
TANZANITE  
BLUE BLOCK

**BLUE BLOCK COATING**



# What is BLUE LIGHT?



- Blue light is a type of visible light with wavelengths between 380-500 nm
- It comes from sunlight, digital screens, and LED lights
- It has high energy and short wavelengths compared to other colors of light

## Blue Light Wavelength Breakdown

↔ **Overall Range**  
380-500 nm

**Harmful Blue-Violet Light**  
400-455 nm


Linked to eye strain and potential eye cell damage, which can affect our circadian rhythm, suppress melatonin production and disrupt our ability to sleep

**Healthy Blue Turquoise Light**  
455-495 nm

- Supports overall well-being by enhancing mood, alertness, and cognitive performance
- Regulates the circadian rhythm, promoting healthy sleep/wake cycles
- Essential for visual performance, aiding color perception and pupil response
- Helps maintain hormonal balance through natural light exposure



**ProctoBLUE Coating is designed to protect your eyes by blocking harmful Blue Violet Light (400 – 455nm) while allowing health supporting Turquoise Light (455 – 495nm) to pass.**



**TINTECH™**


# LENS TINTING











More colour variations  
available under each tint







Expert colour matching  
available upon request





COLOUR VARIATIONS MAY OCCUR DUE TO DIFFERENCE IN LENS MATERIALS



 **INDOOR**



 **CLOUDY**





 **SUNLIGHT**



**ECLIPSAR™** **Sunglare™**

# ECLIPSAR™ Sun glare™




Photochromatic lenses utilize light-sensitive molecules that react to UV exposure, **darkening under sunlight** and **reverting to a clear state indoors or at night**, ensuring optimal visual performance in changing environments.

Also available in SunPink, SunGreen<sup>G15</sup>, SunBlue, and SunBrown, these photochromic lenses offer dynamic light adaptation with a stylish colour finish.


**ECLIPSAR**<sup>TM</sup>  
**Sunpink**

**ECLIPSAR**<sup>TM</sup>  
**Sungreen<sup>G15</sup>**




**ECLIPSAR**<sup>TM</sup>  
**Sunblue**

**ECLIPSAR**<sup>TM</sup>  
**Sunbrown**



# LENS THICKNESS COMPARISON





### Digital Single Vision

<b>82TL021</b>	MYOPIA X WHITE 1.50 UC
<b>82TL025</b>	MYOPIA X SUNGLARE 1.56 HC
<b>82TL023</b>	MYOPIA X WHITE HRI 1.67 HC
<b>82TL016</b>	MYOPIA X SUNGLARE HRI 1.67 HC
<b>82TL033</b>	MYOPIA X WHITE HRI 1.74 HC
<b>83TL001</b>	RELAX ANTI-FATIGUE WHITE 1.50 UC
<b>83TL005</b>	RELAX ANTI-FATIGUE SUNGLARE 1.56 HC
<b>83TL003</b>	RELAX ANTI-FATIGUE WHITE HRI 1.67 HC
<b>83TL006</b>	RELAX ANTI-FATIGUE SUNGLARE HRI 1.67 HC
<b>83TL013</b>	RELAX ANTI-FATIGUE WHITE HRI 1.74 HC
<b>83TL011</b>	EEZ ANTI-FATIGUE WHITE 1.50 UC
<b>83TL015</b>	EEZ ANTI-FATIGUE SUNGLARE 1.56 HC
<b>83TL023</b>	EEZ ANTI-FATIGUE WHITE HRI 1.67 HC
<b>83TL016</b>	EEZ ANTI-FATIGUE SUNGLARE HRI 1.67 HC
<b>83TL033</b>	EEZ ANTI-FATIGUE WHITE HRI 1.74 HC

### Digital Invisible Bifocal

<b>84TL051</b>	TRUE INVISIBLE BIFOCAL FF WHITE 1.50 UC
<b>84TL045</b>	TRUE INVISIBLE BIFOCAL FF SUNGLARE 1.56 HC
<b>84TL023</b>	TRUE INVISIBLE BIFOCAL FF WHITE HRI 1.67 HC
<b>84TL016</b>	TRUE INVISIBLE BIFOCAL FF SUNGLARE HRI 1.67 HC
<b>84TL033</b>	TRUE INVISIBLE BIFOCAL FF WHITE HRI 1.74 HC

### Office Smart

<b>85TL021</b>	OFFICE SMART WHITE 1.50 UC
<b>85TL025</b>	OFFICE SMART SUNGLARE 1.56 HC
<b>85TL033</b>	OFFICE SMART WHITE HRI 1.67 HC
<b>85TL026</b>	OFFICE SMART SUNGLARE HRI 1.67 HC
<b>85TL043</b>	OFFICE SMART WHITE HRI 1.74 HC

### Premium Progressive Lenses

<b>86TL131</b>	ESSENTIAL-Q WHITE 1.50 UC
<b>86TL105</b>	ESSENTIAL-Q SUNGLARE 1.56 HC
<b>86TL163</b>	ESSENTIAL-Q WHITE HRI 1.67 HC
<b>86TL076</b>	ESSENTIAL-Q SUNGLARE HRI 1.67 HC
<b>86TL173</b>	ESSENTIAL-Q WHITE HRI 1.74 HC
<b>86TL141</b>	INFINITY FLEX WHITE 1.50 UC
<b>86TL115</b>	INFINITY FLEX SUNGLARE 1.56 HC
<b>86TL183</b>	INFINITY FLEX WHITE HRI 1.67 HC
<b>86TL086</b>	INFINITY FLEX SUNGLARE HRI 1.67 HC
<b>86TL193</b>	INFINITY FLEX WHITE HRI 1.74 HC
<b>86TL151</b>	MAGIC PRIME 8K WHITE 1.50 UC
<b>86TL125</b>	MAGIC PRIME 8K SUNGLARE 1.56 HC
<b>86TL203</b>	MAGIC PRIME 8K WHITE HRI 1.67 HC
<b>86TL096</b>	MAGIC PRIME 8K SUNGLARE HRI 1.67 HC
<b>86TL213</b>	MAGIC PRIME 8K WHITE HRI 1.74 HC

### Occupational Lens

<b>86TL161</b>	MAX MIDDLE WHITE 1.50 UC
<b>86TL135</b>	MAX MIDDLE SUNGLARE 1.56 HC
<b>86TL223</b>	MAX MIDDLE WHITE HRI 1.67 HC
<b>86TL106</b>	MAX MIDDLE SUNGLARE HRI 1.67 HC
<b>86TL233</b>	MAX MIDDLE WHITE HRI 1.74 HC
<b>86TL171</b>	DRIVE MAX WHITE 1.50 UC
<b>86TL145</b>	DRIVE MAX SUNGLARE 1.56 HC
<b>86TL243</b>	DRIVE MAX WHITE HRI 1.67 HC
<b>86TL116</b>	DRIVE MAX SUNGLARE HRI 1.67 HC
<b>86TL253</b>	DRIVE MAX WHITE HRI 1.74 HC
<b>86TL181</b>	SPORTS AGILITY WHITE 1.50 UC
<b>86TL155</b>	SPORTS AGILITY SUNGLARE 1.56 HC
<b>86TL263</b>	SPORTS AGILITY WHITE HRI 1.67 HC
<b>86TL126</b>	SPORTS AGILITY SUNGLARE HRI 1.67 HC
<b>86TL273</b>	SPORTS AGILITY WHITE HRI 1.74 HC

### Digital Polarised Lenses

<b>86TL027</b>	ESSENTIAL-Q POLARISED 1.50 HC (GREY, BROWN, G15)
<b>86TL037</b>	INFINITY FLEX POLARISED 1.50 HC (GREY, BROWN, G15)
<b>86TL057</b>	MAGIC PRIME 8K POLARISED 1.50 HC (GREY, BROWN, G15)
<b>86TL047</b>	MAX MIDDLE POLARISED 1.50 HC (GREY, BROWN, G15)
<b>86TL067</b>	SPORTS AGILITY POLARISED 1.50 HC (GREY, BROWN, G15)



**PRISMA**  
COLLECTION

**BESPOKE  
PROGRESSIVES**



**OTL LENS LAB**  
YOUR RELIABLE OPTICAL LAB

## CONTACT US

 079 381 9567  [www.prismlensdesigns.com](http://www.prismlensdesigns.com)  
 021 532 1682  [info@prismalensdesigns.com](mailto:info@prismalensdesigns.com)

## NATIONWIDE BRANCHES



**CAPE TOWN**  
27 Moody Ave,  
Epping 1, 7475



**JOHANNESBURG**  
9 Delphi St, Eastgate,  
Sandton, 2090



**DURBAN**  
22 Hillclimb Rd,  
Westmead, 3610



SWISS  
SMART  
INNOVATION

A Product of  
**VISION8**  
GLOBAL



PRISMA™ IS EXCLUSIVELY DISTRIBUTED BY **OTL LENS LAB** IN SOUTHERN AFRICA